



INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(51) International Patent Classification ⁶ : A61L 25/00, A61K 35/16		A1	(11) International Publication Number: WO 98/58689 (43) International Publication Date: 30 December 1998 (30.12.98)
(21) International Application Number: PCT/IT98/00173		(81) Designated States: AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, CA, CH, CN, CU, CZ, DE, DK, EE, ES, FI, GB, GE, GH, GM, GW, HU, ID, IL, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, UA, UG, US, UZ, VN, YU, ZW, ARIPO patent (GH, GM, KE, LS, MW, SD, SZ, UG, ZW), Eurasian patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European patent (AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, ML, MR, NE, SN, TD, TG).	
(22) International Filing Date: 24 June 1998 (24.06.98)			
(30) Priority Data: MI97A001490 24 June 1997 (24.06.97)		IT	
(71)(72) Applicant and Inventor: BERETTA, Roberto [IT/IT]; Via Rho, 8, I-20125 Milano (IT).			
(72) Inventor; and		Published	
(75) Inventor/Applicant (for US only): LODI, Sergio [IT/IT]; Strada per Barbiano, 44, I-22055 Merate (IT).		With international search report. Before the expiration of the time limit for amending the claims and to be republished in the event of the receipt of amendments.	
(74) Agents: PIZZOLI, Pasquale; Società Italiana Brevetti S.p.A., Via Carducci, 8, I-20123 Milano (IT) et al.			

(54) Title: A READY-TO-USE KIT FOR PREPARING AUTOLOGOUS FIBRIN GLUE

(57) Abstract

A ready-to-use kit for preparing autologous fibrin glue comprising a sealed container containing calcium chloride as coagulation activator and, possibly, tranexamic acid or epsilon-amino-caproic acid as fibrin stabilizer.

BEST AVAILABLE COPY

FOR THE PURPOSES OF INFORMATION ONLY

Codes used to identify States party to the PCT on the front pages of pamphlets publishing international applications under the PCT.

AL	Albania	ES	Spain	LS	Lesotho	SI	Slovenia
AM	Armenia	FI	Finland	LT	Lithuania	SK	Slovakia
AT	Austria	FR	France	LU	Luxembourg	SN	Senegal
AU	Australia	GA	Gabon	LV	Latvia	SZ	Swaziland
AZ	Azerbaijan	GB	United Kingdom	MC	Monaco	TD	Chad
BA	Bosnia and Herzegovina	GE	Georgia	MD	Republic of Moldova	TG	Togo
BB	Barbados	GH	Ghana	MG	Madagascar	TJ	Tajikistan
BE	Belgium	GN	Guinea	MK	The former Yugoslav Republic of Macedonia	TM	Turkmenistan
BF	Burkina Faso	GR	Greece	ML	Mali	TR	Turkey
BG	Bulgaria	HU	Hungary	MN	Mongolia	TT	Trinidad and Tobago
BJ	Benin	IE	Ireland	MR	Mauritania	UA	Ukraine
BR	Brazil	IL	Israel	MW	Malawi	UG	Uganda
BY	Belarus	IS	Iceland	MX	Mexico	US	United States of America
CA	Canada	IT	Italy	NE	Niger	UZ	Uzbekistan
CF	Central African Republic	JP	Japan	NL	Netherlands	VN	Viet Nam
CG	Congo	KE	Kenya	NO	Norway	YU	Yugoslavia
CH	Switzerland	KG	Kyrgyzstan	NZ	New Zealand	ZW	Zimbabwe
CI	Côte d'Ivoire	KP	Democratic People's Republic of Korea	PL	Poland		
CM	Cameroon	KR	Republic of Korea	PT	Portugal		
CN	China	KZ	Kazakhstan	RO	Romania		
CU	Cuba	LC	Saint Lucia	RU	Russian Federation		
CZ	Czech Republic	LI	Liechtenstein	SD	Sudan		
DE	Germany	LK	Sri Lanka	SE	Sweden		
DK	Denmark	LR	Liberia	SG	Singapore		
EE	Estonia						

"A READY-TO-USE KIT FOR PREPARING AUTOLOGOUS FIBRIN GLUE"

DESCRIPTION

5 The present invention relates to a kit and a method for preparing autologous fibrin glue, and particularly to a sealed container ready to use being provided with a coagulation activator suitable for obtaining autologous fibrin glue.

10 The fibrin glue is known to be a haemoderivative largely used as a topical surgical adhesive or an haemostatic agent. Several kits are available on the market containing concentrated fibrinogen from donors, associated to a proteic activator of human or animal origin, such as thrombin or batroxobin, for obtaining eterologous fibrin glue.

15 Such known kits involve the use of material of human or animal origin, which, owing to its origin, could result in possible viral contamination and in serious risks for the receiver of the fibrin glue. In the past the authorities have been several times compelled to suspend from the trade or even to ban the haemoderivatives obtained by using material of human or animal origin. Furthermore, rejection cases are known from the literature resulting from the reimplant in the patients of fibrin produced by using human or animal proteins. Such 20 cases are indeed due to the eterologous origin, with respect to the receiver organism, of the sealant protein being reimplanted or some of the components used for preparing it.

25 The autologous fibrin glue, i.e. obtained from the blood of the patient himself, is more reliable with respect to the rejection and/or infection risks. Several procedures have already been described for obtaining extemporary autologous fibrin glue, but no "ready to use" kit is by now available on the market although some relevant references can be found in the patent literature.

30 For example US Patent N.5.733.545 discloses a plasma-buffy coat concentrate to be combined with a fibrinogen activator to form a platelet glue wound sealant. The method disclosed in this patent allows to process blood of the patient for obtaining autologous fibrin glue, but it comprises the use of thrombin or batroxobin as fibrinogen activator. These activators are of human or animal nature and therefore still involve the risk of reject and/or viral infections for the patient.

35 US Patent N. 5.585.007 discloses a method and an apparatus for making concentrated plasma to be used as a tissue sealant. The method consists in

separating plasma from whole blood and removing water from said plasma by contacting it with a concentrator to provide concentrated plasma which can be thereafter coagulated with a solution containing thrombin and calcium. The apparatus comprises a first centrifuge separator in a first chamber, a concentrator 5 included in a second chamber communicating with the first chamber, and a second separator. The method disclosed in this prior art reference requires a long time for obtaining the plasma concentrate necessary for the subsequent preparation of autologous fibrin glue and the apparatus is expensive and not disposable.

The object of the present invention is therefore to provide a ready to use kit, 10 allowing to rapidly obtain autologous fibrin glue and particularly not resulting in viral infections and/or rejection cases when used in surgery.

Such an object is achieved by using a coagulation activator, being neither of human nor of animal origin, but being an inorganic compound which therefore cannot be infected and does not result in rejection.

15 The "ready to use" kit according to the present invention comprises a sealed container containing calcium chloride as coagulation activator. Calcium chloride activates the fibrin present in patient's plasma when this is introduced into the sealed container.

The kit according to the present invention has the great advantage of allowing 20 the preparation of autologous fibrin glue which may be used with no risk of viral infections or rejection cases. Another advantage of the kit according to the present invention is that it allows the preparation of autologous fibrin glue from patient's plasma in a very short time and in the desired form of clots or membrane or spray. Still another advantage of the kit ready to use according to the present invention is 25 to allow the autologous fibrin glue to be obtained at costs proportionally lower with respect to the known systems.

Further advantages of the kit according to the present invention will be evident to those skilled in the art from the following detailed description of some embodiments thereof.

30 As container suitable for the kit according to the present invention can be used for example a glass container for antibiotics as hereinafter described in Example 1. Also glass or plastic test-tubes may be used. The preferred volume of the container is from 5 to 15 ml. The test-tubes have preferably a diameter ranging from 12 to 16 mm and a height ranging from 75 to 100 mm. The container should 35 be suitably thick in order to withstand the stresses resulting from the pressure difference between its inner space and the atmosphere when it is evacuated.

Hemispherical or conical bottom tubes are preferably 0.7 mm thick, flat bottom tubes 1 mm thick. The plastic containers are preferably made of transparent acrylic resin, 0.2-0.8 mm thick, in order to ensure the vacuum keeping for at least 12 months after production. The plastic test-tubes, after the preparation, are preferably 5 introduced into a tin-foil vacuum air-tight container having a heat-sealed inner polyethylene layer, in order to ensure a perfect air-tightness until the date of use.

It should be noted that the evacuation of containers or test-tubes is advisable, however not necessary for putting the present invention into practice.

10 The containers or test-tubes are sealed by rubber or silicon pierceable caps, being suitable to ensure the container to be perfectly air-tight and to allow the vacuum plugging after the introduction of the chemical components and before the steam sterilization step.

15 After the sealing, the containers may be sterilized under steam at 121°C for 30 minutes. The sterilization may be carried out also by irradiation with gamma rays.

As a fibrin stabilizer tranexamic acid can be used, but also pure and crystalline epsilon-amino-caproic acid is suitable to this purpose. The amount will be about 1 g when using a 25 ml container, suitable for a plasma amount of 20 ml. Sometimes it is not necessary to use a fibrin stabilizer.

20 As a coagulation activator solid $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$ is used in the kit according to the present invention. For example, 11.76 mg of $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$ will be introduced in a 5 ml container, by using a precision dosimeter (maximum error: 1-2 mg), in order to prevent polluting foreign components to be introduced.

25 In case of a 15 ml container for a plasma amount of 12 ml, the solid dehydrated calcium chloride amount to be introduced will be as high as 35.28 mg, while the tranexamic acid amount will proportionally be as high as 300 mg of crystals.

30 In case of a 25 ml container for a plasma amount of 20 ml, the dehydrated calcium chloride amount to be introduced will be as high as 58.8 mg, while the tranexamic acid amount will proportionally be as high as 500 mg of crystals.

Besides the dehydrated form used in the Examples, the calcium chloride may be in any other suitable form available on the market, e.g. as $\text{CaCl}_2 \cdot 6\text{H}_2\text{O}$. Also a solution of this salt can be used, as described in Example 1.

EXAMPLE 1

35 In a 5 ml glass container for antibiotics, being sealable under vacuum, made of transparent white glass, inert and 1 mm thick were introduced 100 mg of

tranexamic acid, acting as fibrin stabilizer. The synthetic tranexamic acid, being more than 98% pure, is put on the market by the American company Sigma Inc. Separately, a 1M CaCl_2 solution was prepared, by weighing on a precision balance 147.0 g of $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$ (>99% pure), from the same American company Sigma Inc.

5 This salt was dissolved in exactly 1 liter of ultrapure apyrogenic distilled water, for a few minutes at room temperature, under frequent stirring. By using a precision piston dispenser, having a dispensing precision of $\pm 5\%$ (Eppendorf like), 80 μl of the activator solution were introduced in the glass container. In this step, at the same time as the dispensing, a filtering was carried out by using a 0.22 μm

10 Millipore sterilizing filter, while carefully preventing possible contamination from powders or filaments of any kind. Finally the glass container was plugged with a rubber cap being pierceable and pluggable under vacuum, while minding not to completely plug the container, so as to allow the subsequent vacuum plugging and possibly a further sterilization by using gas. The container was then introduced into

15 a suitable device for vacuum plugging, while preventing any possible contamination from solid particles in the atmosphere (ULPA or HEPA filtration in sterile chamber). A vacuum as high as 4 ml was applied, by using a membrane vacuum pump and a micrometric control, to the inner atmosphere of the device. In order to control the vacuum level in the inner atmosphere, a precision vacuometer was used

20 (precision ± 1 mbar). Finally, without discharging the device, the container was plugged under vacuum, to be thereafter recovered for the use as described in the following Example.

EXAMPLE 2

25 10 ml of venous blood were drawn from a patient according to the provisions of the qualitative standards for clinical analysis, e.g. by using VacutainerTM sterile test-tubes, added with a 0.106 M sodium citrate solution. For this purpose also test-tubes added with disodium or dipotassium etilendiaminetetraacetate can be used. The sample was carefully kept sterile during the blood drawing. Finally, the sample was gently shaken for wholly mixing the components, thereby ensuring the

30 anticoagulating action of sodium citrate. The test-tube was then introduced in a suitable centrifuge, while carefully balancing the rotor weight in order to prevent the same centrifuge to be damaged. Once the lid is sealed, the sample was centrifuged at 3500 rpm for 15 minutes, thereby separating the red cells (being thicker) from the citrated plasma (supernatant). In this case the plasma yield, mainly

35 depending upon the characteristics of the donor blood, was as high as 55%. The test-tube containing the separated plasma was kept plugged in sterile conditions and

was placed vertically in a stand for recovering the plasma itself. In this step care was taken not to shake the test-tube, in order to prevent the mixing of the two phases separated in the centrifugation. The outer portion of the test-tube cap was then sterilized by using denatured alcohol and then a sterile needle, being connected 5 to a sterile syringe, was introduced in the test-tube cap. The needle was brought up to 3-4 mm apart from the separating meniscus of the two phases, and 4 ml of plasma were drawn. By using the same needle, the cap of the container according to the present invention, which had been prepared as described in Example 1, was pierced, having been previously sterilized by using alcohol. As soon as the needle 10 pierced the cap, the citrated plasma contained in the syringe was completely sucked into the container. This was gently shaken and, after about 2 minutes at 37°C, a clot of sterile autologous fibrin glue was obtained, ready to be immediately used.

EXAMPLE 3

About 18 ml of venous blood were drawn from a normotype 49 years-old 15 patient by using 5 ml sodium citrate Vacutainer™ test-tubes, taking care to shake gently just after the drawing of the sample. The so taken blood was immediately subjected to centrifugation (15 min. at 2500 rpm) to separate the plasma. The plasma (12 ml) was carefully transferred into two 10 ml test-tubes, containing 120 µL of CaCl₂ (10g/100ml) each, which had been prepared as described in Example 20 1, but without using tranexamic acid. After mixing the plasma with the activator, the test-tubes were centrifuged for 30 min. at 3000 rpm, finally obtaining two massive fibrin samples which were inserted, with all sterility precautions, within 2-3 hours from preparation, in the large vesicular mandibular cavity resulting from extraction of impacted left canine and right second incisor, as well as from 25 abscission of the cyst present in the central area of the incisor teeth. Finally the gingival edges were closed with eight stitches. A radiographic check 15 days after showed the fibrin still in its position, apparently intact. Histology 7 months after proved the complete replacement of the fibrin with bony tissue, with a better post-operative course than with traditional methods, requiring over 12 months to achieve 30 the same result. Since no antifibrinolytic agent had been used for the preparation of autologous fibrin, it can be stated in this case that said additive was useless for the specific purpose.

EXAMPLE 4

To produce an adhesive fibrin glue 12 ml of plasma, obtained as in Example 35 3, were transferred, with all the measures in order to preserve sterility, into a 20 ml container according to the present invention, prepared as described in Example 1.

After careful stirring, the mixed plasma was poured on a sterile glass slide, of the kind used in chemical laboratories, where the plasma was mixed with sterile and very pure calcium carbonate of coralline origin (BIOCORAL™-INOTEBS S.A. France), or with calcium fluoride (>98% Sigma Inc.). These calcium salts are both 5 well known to the skilled in the art as stimulators of fibroblasts.

By mixing one part of the plasma with one part of calcium carbonate, (e.g.: 2 ml with 500 mg) a malleable, sterile and adhesive paste was obtained and used as a 10 filler for subgingival spaces or different cavities after abscission of infected mucous sacs. The paste, positioned so as to fill the empty spaces, formed in a few minutes a solid fibrin web acting as a haemostatic plug and created an autologous biological substrate supporting the mucous edges in position and where later migration of connectival cells started.

EXAMPLE 5

To obtain a membrane of fibrin glue 20 ml of plasma, obtained as in Example 15 3, were put in a 25 ml, flat-bottomed container according to the present invention prepared as in Example 1. After the usual careful stirring, the container was centrifuged for 40 min. at 4000 rpm with a swing-out rotor. At the end of the centrifuging operation, from the bottom of the test tube a white-colored, very 20 compact and tensile-strong membrane was recovered, having the same size as the bottom of the test-tube (24mm diam.) and thickness of 3 mm. This autologous membrane, owing to its compactness and strength, was used as a holding and separating membrane in dental and general surgery, as a substitute for porous synthetic membranes. The obtained membrane can be sterile stored for several days at 4°C.

EXAMPLE 6

To obtain large-sized membranes of fibrin glue about 200 ml of citrated plasma were drawn from a patient, collected and separated in a double transfusion bag. The plasma was subjected to cryoprecipitation by freezing at -80°C for 12 30 hours, defreezing being carried out overnight at 4°C (this procedure is well known to those skilled in the art). The same morning the plasma obtained by this procedure was subjected to centrifugation for 15 min. at 5000 rpm at 4°C to obtain about 20 ml of cryoprecipitate. After careful removal of the supernatant by using a pressing device (e.g. XP100 of the company Jouan S.A. France) the cryoprecipitate was taken up with 20 ml of whole plasma of the same patient. The resulting 40 ml were 35 put in a 35 mm diameter, flat-bottomed sterile polypropylene container according to the present invention, containing the suitable quantity of activator, as described

in Example 1. After careful shaking, the container was centrifuged for 40 min. at 5000 rpm to obtain a membrane as in Example 5, but more compact and tensile-strong owing to the higher content of fibrin. Said membrane too can be stored in sterile form for several days at 4°C.

5 The membrane obtained by the method described in Example 5, in addition to utilization described in Example 4, can be used as a substrate for the culture in vitro of dermal cells of the same patient, in order to obtain grafts to be transplanted in case of very serious scalds.

10 Membranes of a good quality useful for the above mentioned purposes can be obtained also from whole separated plasma directly transferred into the container according to the present invention. The obtained membrane will be thinner than the above described one, but still useful for surgical uses and as a substrate for cellular growth.

EXAMPLE 7

15 To obtain spray fibrin starting from a cryoprecipitate as in Example 5, 20 ml of cryoprecipitate were taken up with 10 ml of whole plasma at room temperature and gently shaken, to complete dissolution. The resulting plasma was carefully transferred into a 50 ml container according to the present invention prepared as in Example 1, shaking gently for a perfect mixing of the components. After 120 sec. at 20 room temperature, the test-tube was connected to a Venturi-type sterile air compressor, known to those skilled in the art, to be uniformly distributed on the surface of a bleeding organ being subjected to surgery (lung, heart, spleen, arterious anastomosis). The concentrated plasma, containing concentrated fibrinogen, thrombin, calcium ions and other coagulation enzymes, distributed over the organ, 25 coagulated within a few seconds, owing also to tissue coagulation activating enzymes present in the endothelium of the patient creating a fibrin film having a protective haemostatic function. The surgical operation was therefore concluded with the reduction of internal haemorrhages and so avoiding further blood transfusions or complications.

CLAIMS

1. A ready to use kit for preparing autologous fibrin glue characterized in that it comprises a sealed container containing calcium chloride as coagulation activator.
- 5 2. A kit according to claim 1, characterized in that the sealed container further contains a fibrin stabilizer.
3. A kit according to claim 2, characterized in that tranexamic acid is used as fibrin stabilizer.
- 10 4. A kit according to claim 2, characterized in that epsilon-amino-caproic acid is used as fibrin stabilizer.
5. A kit according to any of the previous claims, characterized in that the sealed container is under vacuum.
- 15 6. A method for preparing autologous fibrin glue by means of the kit according to any of the previous claims, characterized in that it comprises drawing the blood from a patient, separating plasma from blood and contacting plasma with calcium chloride in a sealed container under a sterile atmosphere.

1
INTERNATIONAL SEARCH REPORT

International application No.

PCT/IT 98/00173

A. CLASSIFICATION OF SUBJECT MATTER

IPC6: A61L 25/00, A61K 35/16

According to International Patent Classification (IPC) or to both national classification and IPC

B. FIELDS SEARCHED

Minimum documentation searched (classification system followed by classification symbols)

IPC6: A61L, A61K

Documentation searched other than minimum documentation to the extent that such documents are included in the fields searched

Electronic data base consulted during the international search (name of data base and, where practicable, search terms used)

WPI, EPODOC, CA, MEDLINE, BIOSIS, EMBASE

C. DOCUMENTS CONSIDERED TO BE RELEVANT

Category*	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
X	WO 9627397 A1 (QUANTIC BIOMEDICAL PARTNERS), 12 Sept 1996 (12.09.96), page 8, line 16 - page 9, line 6, the abstract --	1-6
X	WO 9617871 A1 (PLASMASEAL CORPORATION), 13 June 1996 (13.06.96), page 33, line 15 - line 20; page 34, line 12 - line 17 --	1-6
X	WO 9512371 A1 (CLARION PHARMACEUTICALS, INC.), 11 May 1995 (11.05.95), page 3, line 2 - line 10, the abstract --	1-6

 Further documents are listed in the continuation of Box C. See patent family annex.

- * Special categories of cited documents:
- "A" document defining the general state of the art which is not considered to be of particular relevance
- "E" earlier document but published on or after the international filing date
- "L" document which may throw doubts on priority claim(s) or which is cited to establish the publication date of another citation or other special reason (as specified)
- "O" document referring to an oral disclosure, use, exhibition or other means
- "P" document published prior to the international filing date but later than the priority date claimed
- "T" later document published after the international filing date or priority date and not in conflict with the application but cited to understand the principle or theory underlying the invention
- "X" document of particular relevance: the claimed invention cannot be considered novel or cannot be considered to involve an inventive step when the document is taken alone
- "Y" document of particular relevance: the claimed invention cannot be considered to involve an inventive step when the document is combined with one or more other such documents, such combination being obvious to a person skilled in the art
- "&" document member of the same patent family

Date of the actual completion of the international search

24 Sept 1998

Date of mailing of the international search report

28.10.98

Name and mailing address of the ISA/
European Patent Office, P.B. 5818 Patentlaan 2
NL-2280 HV Rijswijk

Authorized officer

GADOL TMA DAI McDANTZ

INTERNATIONAL SEARCH REPORT

International application No.

PCT/IT 98/00173

C (Continuation). DOCUMENTS CONSIDERED TO BE RELEVANT

Category*	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
X	EP 0592242 A1 (E.R. SQUIBB & SONS, INC.), 13 April 1994 (13.04.94), page 4, line 27 - line 33; page 5, line 1 - line 2; claim 64 --	1-6
X	Arch Otorhinolaryngol, Volume 237, 1983, G. Wolf, "Der konzentrierte autologe Gewebekleber", page 279 - page 283; page 280; paragraphs 3-6; page 281, paragraphs 1-3 --	1-6
A	US 5030215 A (BRENDA S. MORSE ET AL), 9 July 1991 (09.07.91), column 2, lines 2-7, 43-45; column 3, lines 48-51; and column 4, lines 42-48, 59-60 --	1-6
A	WO 9422503 A1 (OPPERBAS HOLDING B.V.), 13 October 1994 (13.10.94), the claims -- -----	3

INTERNATIONAL SEARCH REPORT

Information on patent family members

27/07/98

International application No.

PCT/IT 98/00173

Patent document cited in search report	Publication date	Patent family member(s)	Publication date
WO 9627397 A1	12/09/96	AU 5416696 A CA 2214274 A EP 0813427 A US 5733545 A	23/09/96 12/09/96 29/12/97 31/03/98
WO 9617871 A1	13/06/96	AU 4286396 A BR 9506744 A CA 2182862 A EP 0743957 A JP 9509432 T NZ 297652 A US 5585007 A US 5788662 A	26/06/96 16/09/97 13/06/96 27/11/96 22/09/97 26/05/97 17/12/96 04/08/98
WO 9512371 A1	11/05/95	AU 1086795 A CA 2175203 A EP 0726749 A JP 9504719 T US 5643596 A US 5645849 A	23/05/95 11/05/95 21/08/96 13/05/97 01/07/97 08/07/97
EP 0592242 A1	13/04/94	AU 676201 B AU 4887893 A BR 9304185 A CA 2108104 A CN 1091315 A CZ 9302117 A ES 2065294 A,B HU 66587 A HU 9302839 D IL 107211 A JP 6199685 A MX 9306272 A NO 933624 A NZ 248897 A NZ 280264 A PL 173858 B SK 109393 A US 5739288 A US 5750657 A US 5763410 A US 5763411 A US 5770194 A US 5773418 A	06/03/97 21/04/94 03/05/94 09/04/94 31/08/94 17/04/96 01/02/95 28/12/94 00/00/00 30/09/97 19/07/94 29/07/94 11/04/94 27/08/96 22/08/97 29/05/98 08/06/94 14/04/98 12/05/98 09/06/98 09/06/98 23/06/98 30/06/98
US 5030215 A	09/07/91	AU 7220691 A CA 2072371 A EP 0509054 A JP 5504280 T WO 9109573 A	24/07/91 04/07/91 21/10/92 08/07/93 11/07/91

INTERNATIONAL SEARCH REPORT

Information on patent family members

27/07/98

International application No.

PCT/IT 98/00173

Patent document cited in search report	Publication date	Patent family member(s)	Publication date
WO 9422503 A1	13/10/94	AT 151295 T	15/04/97
		AU 4566893 A	14/02/94
		AU 6536394 A	24/10/94
		CA 2140427 A	03/02/94
		CA 2159469 A	13/10/94
		CZ 9500054 A	14/06/95
		DE 69309685 D,T	16/10/97
		EP 0651659 A,B	10/05/95
		EP 0691858 A	17/01/96
		FI 950194 A	17/01/95
		HU 71873 A	28/02/96
		HU 9500141 D	00/00/00
		IL 109157 D	00/00/00
		JP 8502035 T	05/03/96
		JP 8512215 T	24/12/96
		NO 950171 A	17/01/95
		SK 5095 A	11/07/95
		WO 9402182 A	03/02/94

**This Page is Inserted by IFW Indexing and Scanning
Operations and is not part of the Official Record**

BEST AVAILABLE IMAGES

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

- BLACK BORDERS**
- IMAGE CUT OFF AT TOP, BOTTOM OR SIDES**
- FADED TEXT OR DRAWING**
- BLURRED OR ILLEGIBLE TEXT OR DRAWING**
- SKEWED/SLANTED IMAGES**
- COLOR OR BLACK AND WHITE PHOTOGRAPHS**
- GRAY SCALE DOCUMENTS**
- LINES OR MARKS ON ORIGINAL DOCUMENT**
- REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY**
- OTHER:** _____

IMAGES ARE BEST AVAILABLE COPY.

As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.